

BlackBox™

Manual Endoscope Processing System

A bright new concept from STERIS




BlackBox delivers to
manual endoscope
processing:

- > Ease of use
- > Direction
- > Traceability
- > Rigour
- > Repeatability


Manual processing of endoscopes is an important part of hospital infection control throughout the world.

BlackBox is designed to help the busy hospital perform the right processing steps every time. BlackBox uses icons to make the process easy, and a printer to record each step for traceability. The endoscope processing steps include:



- > The hospital selects which steps to follow.
- > BlackBox guides the operator with icons.
- > The printer records the process.

STERIS has been a leader in Infection Prevention for over 100 years. Careful thought and generations of practical experience have created BlackBox: a modern, smart, easy to use solution for better practice



STERIS Corporation
5960 Heislley Road
Mentor, OH 44060-1834 ■ USA
440-354-2600 ■ 800-548-4873
www.steris.com

STERIS Limited
STERIS House
Jays Close
Viabiles
Basingstoke
Hampshire RG22 4AX ■ UK
+44 1256 840400 ■ 0800 252609
Email: UKHCsales@steris.com
www.steris.com/international



STERIS OFFICES WORLDWIDE	
Brazil	81 3 5521 2104
Canada	32 3 389 1197
China	1 805 677 0863
France	33 2 38 70 83 50
Germany	49 221 466120 30
Greece	30 210 6800 848
India	91 33 2387 5150
Italy	39 022 13 034 1
Japan	81 3 5521 2104
Korea	82 2 517 1517
Latin America	1 800 884 9550
Malaysia	6 03 7954 9822
Singapore	65 68 41 7677
Spain	34 91 658 5920
Switzerland	41 32 376 0200
United Kingdom	44 1256 840400

Document # M3470GB 2009.11, Rev. A
©2011 STERIS Corporation.
All rights reserved.
TPS Printed 03/2011, 1000

Technology to Prevent Infection and Contamination™